


Wall System AW421




AW4211 AW4213



AW4212 with left stop AW4212 with right stop



Base plate



Base plate

- AW4211** Wall system for 1 article of clothing
- AW4213** Wall system for 3 articles of clothing

- AW4212R** Wall system for 2 articles of clothing *with right stop*
- AW4212L** Wall system for 2 articles of clothing *with left stop*

AW421 – Storage System for X-Ray Protective Clothing

The wall system with top marks for user-friendliness and storage space.

Up to three swivelling chrome bars offer you a range of operation of up to 180°. In order to create space in the room, a maximum of two bars can be folded to one side. The systems AW4212R (right) and AW4212L (left) differ in the side of the stop position. In the case of the system AW4213 for three aprons, the bars have to be spread to the right and the left side.

The base plate is painted in medico-white (RAL 9002), the swivel bars are chromeplated.

Technical Data

	AW4211	AW4212R	AW4212L	AW4213
Weight	1.2 kg	1.7 kg	1.7 kg	2.2 kg
Max. Load	12 kg	24 kg	24 kg	36 kg
Hanging Capacity	1*	2*	2*	3*
Dimensions (base plate)	205 x 85 mm	205 x 85 mm	205 x 85 mm	205 x 85 mm
Projection (swivel bar)	630 mm	630 mm	630 mm	630 mm
Colour	medico-white (RAL 9002)	medico-white (RAL 9002)	medico-white (RAL 9002)	medico-white (RAL 9002)

*X-ray protective aprons / coats / suits

MAVIG GmbH
X-Ray Protection
& System Solutions

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone: +49 (0) 89 / 420 96 – 0
Fax: +49 (0) 89 / 420 96 – 200
e-Mail: info@mavig.com

www.mavig.com