

# CM TraumaCover™

CM ProLine



Fast and stable fixation of the head.



Every package comes with a complete head fixation.

- ✓ WARMING
- ✓ HYGIENIC
- ✓ TIMESAVING
- ✓ SAFE  
IMMOBILISATION

#### PRODUCT INFORMATION:

CM TraumaCover is a single use cover for traumamattresses.

The material is enviromentally friendly, consisting of polyurethane and polyester. It is waterresistant with a high tensile strength.

The head can be easily and quickly immobilized using the disposable head cushion.

Two velcro straps for fixation of the head and a head cushion come with each CM Trauma Cover.

Size: 225 x 162 cm.

- Art.Nr: CMCOVER
- Name: CM TraumaCover
- Package Quantity: 15 pcs

Made in Sweden by Comfort Medical AB.



DISTRIBUTOR: