

# CT-CUSHIONS

CM ProLine



Hard / soft cushion 40 cm.


Welded covers can be wiped with alcohol and soap solution.

Also available in 25 cm length, suitable for knee pads.

Article nr: CMPLKH40 (hard, 40x17x8 cm)

CMPLKM40 (soft, 40x17x8 cm)

CMPLKM25 (soft, 25x17x8 cm)



Wedge for CT headrest.

Available in two heights, 3 cm and 5 cm.

Made in EVA, can be wiped with alcohol / soap solution.

Article nr: CMPCTKIL3 ( 3 cm height)

CMPCTKIL5 ( 5 cm height)




Pressure relieving cushion for CT headrest.

PU-coating, can be wiped with alcohol / soap solution.

size: 30x18 cm

Art.nr: CMPCTKUDDE




Tiltable PU-coated cushion for headrest.

Suitable for examination outside of the headrest.

Can be wiped with alcohol / soap solution.

Articlenr: CMPHD



OverHead CT is used for positioning the patients arms at CT examinations. Can also be used as a knee cushion.

Dimension: 43x30 cm, height 26 cm. Can be wiped with alcohol / soap solution.

Art.nr: OVERH



Kneecushions

Coating: Technoskin, artificial leather.

Core: Polyether

CMPKKF flat 55x38, height 16-8 cm

CMPKK 45x31 cm, height 17 cm

CMPKKC Cylinder 55x15 cm

**DISTRUBUTOR:**

**COMFORT** medical

**Manufactured by:**

Comfort Medical AB

[www.comfortmedical.se](http://www.comfortmedical.se)