

RP278 – Ovary Shield

Extensive tests created the basis for the size and positional relationship between the fixation point "Spinae iliacae anteriores" and the ovarian protection applied. The two radiolucent centring marks ensure correct and simple fixation.

A complete set (all 3 sizes) enables the best possible protection for all your patients. In addition, each model can be modified by means of a "fan" principle, ensuring optimal fine adjustment.

The RP278 is, of course, easy to clean and disinfect.



The design of this MAVIG ovary shield is unsurpassed and meets all requirements of radiology practice. Radiation exposure of the female gonads is considerably reduced without impairing the diagnostic procedure.

RP276 – Testicle Shield for Boys/Men

This model for the protection of male gonads also completely covers the penis and scrotum with a lead equivalent value of Pb 1.00 mm.

MAVIG has, of course, made sure that using the shield does not cause any risk of injury to the patient.

The hook & burr closure provides good stability when adjusted correctly. This model can be unfolded by means of snap fasteners for cleaning and disinfection.

Our range also includes matching disposable hygiene bags.

RP280 – Testicle Shield for Boys/Men

The penis and scrotum are completely protected with a lead equivalent value of Pb 1.00 mm.

A padded, elastic case clasp ensures a safe and correct fit with gentle pressure for the duration of the entire test. Due to the rounded design there is no risk of injury to the patient.

Putting it on is simple and the skin-friendly coating of the MAVIG RP280 testicle capsules provides perfect hygiene inside and out. Cleaning and disinfection are easy.

Of course, radiation protection is the focus of our product design and with the use of the RP280 series, you can guarantee your patients that perfect protection of the genital area in compliance with all legal specifications is ensured.



Protective Material

Standard lead
(RP278 and RP276)

Single layer lead rubber
(RP280)

Lead Equivalent

Standard: Pb 1.00 mm

Outer Material / Colours

Soft PVC cover in light grey
(RP278 and RP276)

Skin-friendly coating in medico-white
(RP280)

Sizes

RP2781 Large height 12 cm
RP2782 Medium height 9 cm
RP2783 Small height 7 cm

RP276M For men, dimensions approx. 155 x 120 x 90 mm
RP276K For boys, dimensions approx. 110 x 90 x 65 mm

RP280M For men, dimensions approx. 135 x 105 x 60 mm
RP280K For boys, dimensions approx. 105 x 85 x 55 mm

Price Saving

The ovary shield is also available as a complete set with all 3 sizes (RP278SET).

Accessories

Hygienic Bag for RP276/RP280

- RP-HYG-CM - Set of 100 bags in men's size
- RP-HYG-CK - Set of 100 bags in boys' size
- RP-HYG-CMK - Set of 100 bags in men's size and 100 bags in boys' size

Our patient protection products conform to the medical devices directive 93/42/EEC.

MAVIG X-ray protective clothing gives you the security of fulfilling all applicable legal requirements.