

CM X-Ray Mattress

CM ProLine



We deliver mattresses to several leading x-ray companies.

PRODUCT INFORMATION:

CM X-ray Mattress is a comfortable mattress with eight handles. The mattress can be custom made in different sizes.

The core of the mattress consists of viscoelastic foam and EVA-foam. The cover is made of durable medical PU-fabric. The mattress underside have low friction attributes for easy adjustments when needed. The standard thickness is 3 cm.

- **COMFORTABLE**
- **DURABLE**
- **ERGONOMIC**

Article no: CM + width and length
For example: CM65220

Cleaning and disinfection:

Disinfect with alcohol and/or soap solution. Remove blood stains as soon as possible with cold water.

Safety:

CM X-ray mattress has an underside with low friction, be aware when transferring patients from and to the mattress.

CE

Last updated 2020-07-23