



## Dental Apron for Images When Sitting

- **Size:** Universal
- **Lead Equivalent:** Pb 0.50 mm
- **Inner Material:** Multi-Layer Lead Rubber
- **Outer Cover:** ComforTex® HPMF  
Dual Coloured & Reversible

**RD64412** Dental Apron with Collar, Dual Coloured: Curacao and Orchidee  
**RD64457** Dental Apron with Collar, Dual Coloured: Limette and Ocean  
 Dimensions: Length 90.0 cm; Width 60.0 cm (not including collar)



## MAVIG X-Ray Protective Dental Apron with Collar

MAVIG's products offer solid protection and optimized functionality. Our new dental apron comes with a specially designed collar for a comfortable fit combined with optimized protection of the thyroid and a quick & easy closure thanks to a velcro fastener on the neck.

The generous cut of this dental apron - for images when sitting - protects the whole shoulder area and two added loops on the back can be attached to a dental chair. Extra features are the universal size and that it is easy to be put on and take off.

### Inner Material: Multi-Layer Lead Rubber

Our natural lead rubber - without any artificial plasticisers - has kept us in the lead as manufacturers of X-ray protective clothing for many years, with excellent marks for weight distribution and radiation protection.

### Outer Cover: ComforTex® HPMF - Dual Coloured & Reversible

Our ComforTex® HPMF, a high performance medical fabric according to DIN EN 13795, provides a far greater protection of the radiation protective inlay than ever before. Its improved barrier function complies with the very demanding hygienic standards in today's work environment.

### MAVIG GmbH

X-Ray Protection  
& System Solutions

Stahlgruberring 5  
81829 Munich  
Germany

PO Box 82 03 62  
81803 Munich  
Germany

Phone: +49 (0) 89 420 96 - 0  
 Fax: +49 (0) 89 420 96 - 200  
 e-Mail: info@mavig.com

[www.mavig.com](http://www.mavig.com)