


# MRI CUSHIONS

CM ProLine



- ✓ FORMBABLE CUSHIONS
- ✓ HYGIENIC
- ✓ EASY TO CLEAN

MRI cushions for head stabilization. This set with four cushions give the patient good support for the head. The cushions is moldable and provides a certain hearing protection.



## PRODUCT INFORMATION

Article nr: CMPMRO

Coating: PU material

Core: Viscoelastic material

Size: 14x9x7/2 cm

Article nr: CMPMRHS

Coating: PU material

Core: Viscoelastic material

Size: 14x9x2 cm

Cleaning: Disinfection with alcohol/soap solution.

**DISTRUBUTOR:**

**COMFORT** medical

**Manufactured by:**  
Comfort Medical AB  
[www.comfortmedical.se](http://www.comfortmedical.se)