

SANDBAGS



Article no: SA0103
Size: 35x13 cm
Weight: 2 kg



Article no: SA0102
Size: 70x10 cm
Weight: 2,5 kg



Artikel nr: ODE18-03
Size: 40x15 cm, Height 5 cm
Weight: 2,3 kg
100% MRI safe positioning bag.
Adapted for use in magnetic camera examination.

- **FLEXIBLE USE**
- **HYGIENIC**

PRODUCT DESCRIPTION

Sandbags in a varied range can give support and facilitate examinations. The bags are equipped with double covers. The outer cover can easily be disinfected.

The design of the bags offers flexibility and the opportunity to be placed as needed.

Cleaning instructions:

Disinfect with alcohol and / or soap solution. Remove blood stains as soon as possible with cold water.

Manufactured by MID.

Last updated 2020-10-06