

„State of the Art Eye Protection“



BR126 - Design: Cocoa / Light Blue

BR126 - Design: Silver Bronze / Orange

Special features of the BR126 & BR127 models:

- Great frame stability through high-grade synthetic materials
- Optimized balance of sturdy, light-weight frame and lenses
- No plasticizers used throughout the entire manufacturing process and low allergy risk due to carefully chosen materials
- Minimized risk of injury due to shape and used components
- Easy cleaning and superior resistance against chemical influences (disinfection)
- Robust and durable finish coating

BR126

The model BR126 offers the choice of two frame sizes: *small-medium or large*. The design of the glasses adapts harmoniously to the facial anatomy resulting in a comfortable fit. The model design features extended lateral protection through additional inlays in the frame.

Included in delivery: case, headband, and lens cleaning cloth

**Lead equivalent values:** *Front/side protection*  
 (weight) Pb 0.50 mm / 0.50 mm (approx. 70 g)  
 Pb 0.75 mm / 0.50 mm (approx. 80 g)

This model comes standard with antireflective coated lenses.

This model is also available with *single or progressive lenses*.

MAVIG only uses protective lenses of the highest quality. This guarantees the best light transmission values. An ideal composition of the protective material results in a high level X-ray absorption and optimized production processes achieve an absolute homogeneous allocation of the glass additives.



Dimensions:	Width of Frame	Lens Height / Lens Width	Bridge Width	Temple Length
	140 mm (S/M)	36 mm / 54 mm (S/M)	17 mm (S/M)	130 mm (S/M)
	145 mm (L)	36 mm / 58 mm (L)	16 mm (L)	130 mm (L)