

# CM T-Fix™

CM ProLine



Fast and stable fixation of the head.



Every CMTFIX comes with a complete head fixation.

- ✓ HYGIENIC
- ✓ TIMESAVING
- ✓ SAFE IMMOBILISATION

#### PRODUCT INFORMATION:

CM T-Fix is a single use head fixation for trauma mattresses. Every package comes with one cushion, two straps for the head and a girdle for fixation. The material is environmentally friendly, consisting of polyurethane and polyester.

The girdle is waterresistant.

The head can be easily and quickly immobilized using the disposable polyester head cushion. The cushion is attached with velcro direct to the girdle.

Art.nr: CMTFIX

Name: CM T-Fix Head Fixation

Package Quantity: 25 pcs

Made in Sweden by Comfort Medical AB



Single use only

DISTRIBUTOR: