

CM BACK CUSHION

CM ProLine



PRODUCT DESCRIPTION

The back cushion can be used to prevent the patient from sinking into the mattress while doing X-ray pulm side. Two handles facilitate handling.

The cushion is stable with a top layer of EVA foam. Below is a core of polyether. This design avoid the patient from sinking down into the cushion.

- **HYGIENIC**
- **COMFORTABLE**
- **STABLE**

Article no: CMPRK

Size: 55x36 cm, height 16-8 cm

Material:

Cover: Artificial leather

Core: Polyether foam / EVA foam

Cleaning

instructions: Disinfect with alcohol and / or soap solution. Remove blood stains as soon as possible with cold water.

Safety instructions:

- The cushion has limited flexibility and should be removed as soon as possible after examination.



Last updated 2020-10-05