

CM ArmFix

CM ProLine



Article no: CMPARM

- **STABLE**
- **EASY TO USE**
- **COMFORTABLE**

Article no: CMPARM

Size: 45x20 cm, Height 37 cm

Cover: Artificial leather (PU)

Core: Polyether foam

Cleaning instructions:

Disinfect with alcohol and / or soap solution. Remove blood stains as soon as possible with cold water.

CE

2021-04-28 (Edition:2)

Distributor: **COMFORT** medical

Manufacturer: Bima Plastteknik AB,
Trålvägen 5, 43274, Träslövsläge, Sweden

Comfort Medical AB | Hjalmar Brantings väg 11G | 26143
Landskrona | Sweden | +46 (0)733-109020 |
Info@comfortmedical.se | www.comfortmedical.se