

# CM X-Ray Mattress

CM ProLine



We deliver mattresses to several leading X-ray companies.

## PRODUCT INFORMATION

CM X-Ray Mattress is a comfortable mattress manufactured in different sizes and thicknesses in order to fulfill a varied demand of X-ray mattresses.

CM X-Ray mattress is manufactured in two different thicknesses (2,5 and 3 cm). The mattresses with a thickness of 3 cm are manufactured with a surface of PU-fabric, a underside with low friction attributes and eight handles which facilitates simple adjustments when needed.

The Mattresses with a thickness of 2,5 cm are mainly intended for use in MRI-scanner. These mattresses are manufactured with PU-fabric without low friction underside or handles.

All mattress cores consists of a combination of viscoelastic polyurethane foam and a layer of EVA-foam.

● **COMFORTABLE**

● **DURABLE**

● **ERGONOMIC**

### Article numbers

See next page

### Cleaning and disinfection

Disinfect with alcohol and/or soap solution. Wipe clean the mattress with water after use of disinfectants. Remove blood stains as soon as possible with cold water.

### Safety instructions

- Be careful when transferring patients to and from the X-ray mattress.
- CM X-Ray mattress are not made for lifting patients on the mattress.
- Consider patient individual risk of developing pressure ulcers in case of long examination times.
- Clean and disinfect the mattress according to given instructions in between patients.
- Avoid contact with warm sources or sharp objects, it can damage the mattress.
- Examine the mattress regularly for damage. Damaged mattress should not be in use.

CE

2021-03-15 (edition:4)

Distributor: **COMFORT** medical

Comfort Medical AB | Hjalmar Brantings väg 11G | 26143  
Landskrona | Sweden | +46 (0)733-109020 |  
Info@comfortmedical.se | www.comfortmedical.se

Manufacturer: Bima Plastteknik AB, Trålvägen 5, 43274, Träslövsläge, Sweden

# CM X-Ray Mattress

CM ProLine

## ARTICLE NUMBERS

- **X-Ray Mattresses with a thickness of 3 cm, top surface of PU-fabric with low friction underside and eight handles.**

Article number: CM+Width+Length

\*Exception CMADORA

Length/Width	57 cm	60 cm	65 cm	68,5 cm	70 cm	75 cm
200 cm	CM57200	CM60200	CM65200		CM70200	CM75200
208,5 cm				CMADORA		
210 cm	CM57210	CM60210	CM65210		CM70210	CM75210
220 cm	CM57220	CM60220	CM65220		CM70220	CM75220
230 cm	CM57230	CM60230	CM65230		CM70230	CM75230
240 cm	CM57240	CM60240	CM65240		CM70240	CM75240

- **MRI Mattresses with a thickness of 2,5 cm, manufactured with PU-fabric, without low friction underside and no handles.**

Article number: CMP+Width+Length+MR

Length/Width	50 cm	52 cm	55 cm
155 cm	CMP50155MR	CMP52155MR	CMP55155MR
165 cm	CMP50165MR	CMP52165MR	CMP55165MR
175 cm	CMP50175MR	CMP52175MR	CMP55175MR
180 cm	CMP50180MR	CMP52180MR	CMP55180MR