


PLASTER CUSHIONS

CM ProLine



Article no: CMPRFEMT

Article no: CMPRFEMT25

PRODUKTBESKRIVNING

Plaster cushions are manufactured with a core of massive sensaten. The cover is made of artificial leather.

The cushions have good stability and can be used for plastering the lower legs.

- **STABLE**
- **EASY TO CLEAN**

Article no: CMPRFEMT25

Size: 25x18x17 cm

Cover: Artificial leather

Core: Sensaten foam

Article no: CMPRFEMT

Size: 38x18x17 cm

Cover: Artificial leather

Core: Sensaten foam

Cleaning

instructions: Disinfect with alcohol and / or soap solution. Remove blood stains as soon as possible with cold water.

Recycling: No disassembly required. Sort as combustible waste.



Last updated 2020-10-05