

# Positioning Cushions

CM ProLine



## PRODUCT INFORMATION

The larger wedge-shaped positioning cushions are made of massive polyether with a layer of EVA and a PU-cover. The cushions are stable and serve, for example, as a good support when positioning the patient on the side.

The smaller cushions are mainly designed as a general support/positioning tools during X-ray examinations.

- **X-RAY COMPATIBLE**

- **EASY TO CLEAN**

### Article no

For more information see next page.

### Cleaning instructions

Disinfect with alcohol and / or soap solution. Remove blood stains as soon as possible with cold water.

### Recycling

No disassembly required. Sort as combustible waste.

Last updated 2022-04-06 (Edition:3)


Distributor: **COMFORT** medical

Manufacturer: Bima Plastteknik AB,  
Trålvägen 5, 43274, Träslövsläge, Sweden


Comfort Medical AB | Hjalmar Brantings väg 11G | 26143  
Landskrona | Sweden | +46 (0)733-109020 |  
Info@comfortmedical.se | www.comfortmedical.se

# Positioning Cushions

CM ProLine




**Article no:** CMP60S  
Wedge-shaped cushion  
60 x 25 x 13 cm, 45°,  
with PU-cover




**Article no:** CMP50S  
Wedge-shaped cushion  
50 x 25 x 13 cm, 45°,  
with PU-cover




**Article no:** CMP40S  
Wedge-shaped cushion  
40 x 25 x 13 cm, 45°,  
with PU-cover




**Article no:** CMPF40S  
Wedge-shaped cushion  
40 x 24 x 12 cm, 30°,  
with PU-cover




**Article no:** CMPPB  
Foot cushion for  
follow-up after hip surgery  
24/12 x 18 x 25 cm




**Article no:** CMP2040  
Track cushion  
40 x 20 x 5 cm




**Article no:** CMPHK  
Hand cushion  
18 x 15 x 6 cm



**Article no:** CMPHT  
Stair cushion for fingers  
13 x 13 x 8 cm



**Article no:** CMPKS  
Low wedge  
19 x 19 x 0,5-4 cm 10°



**Article no:** CMPKL  
High wedge  
19 x 19 x 0,5-6 cm 16°



**Article no:** CMPHLK/CMPHLKS  
Wrist wedges  
23x16x10 / 18x9x4,5 cm 45°