

Knee Cushions

CM ProLine



Article no: CMPKK
Name: Knee cushion
Size: 45x31cm
Height: 17 cm
Cover: PU-cover
Core: Polyether foam



Article no: CMPKKF
Name: Knee cushion flack
Size: 55x38 cm
Height: 16-8 cm
Cover: PU-cover
Core: Polyether foam



Article no: CMPKKC
Name: Cylinder knee cushion
Size: 50x15 cm
Cover: PU-cover
Core: Polyether foam



Article no: CMPLKM25
Name: Small knee cushion 25
Size: 25x17cm
Height: 8 cm
Cover: PU-cover
Core: Polyether foam

Article no: CMPLKM40
Name: Small knee cushion 40
Size: 40x17 cm
Height: 8 cm
Cover: PU-cover
Core: Polyether foam

Article no: CMPLKH40
Name: Small knee cushion 40 (Hard)
Size: 40x17 cm
Height: 8 cm
Cover: PU-cover
Core: Sensaten

● **COMFORTABLE**

● **ERGONOMIC**

PRODUCT DESCRIPTION

CM Proline offer a variety of comfortable knee cushions. Main objectives and purpose are hygiene, ergonomics and comfort. Knee cushions are available in a range of shapes to fullfil different needs of positioning. All knee cushions are covered in outer material that can easily be cleaned.

CE

2022-04-01 (Editon:2)

Distributor:

COMFORT medical

Manufacturer: Bima Plastteknik AB, Trålvägen
5, 43274, Träslövsläge, Sweden

Comfort Medical AB | Hjalmar Brantings väg 11G | 26143
Landskrona | Sweden | +46 (0)733-109020 |
Info@comfortmedical.se | www.comfortmedical.se