

CM OverHead CT

CM ProLine



PRODUCT DESCRIPTION

CM Overhead CT is a cushion made for fixation of arms while doing CT thorax/abdomen. The cushion provide a stable and restful position for the patients arms. The cushion is made with a core of polyether foam and a PU-cover. The cushion comes with two straps which can be strategically placed to fixate the patients arms.

- STABLE FIXATION
- EASY TO USE
- COMFORTABLE

Article number

OVERH

Size

43x30 cm, Height 26 cm

Material

Core: polyether foam

Cover: PU-Cover

CE

2022-04-01 (Utg:2)

Distributör:

COMFORT medical

Manufacturer: Bima plastteknik AB, Trålvägen 5,
43274, Träslövsläge, Sweden

Comfort Medical AB | Hjalmar Brantings väg 11G | 26143
Landskrona | Sweden | +46 (0)733-109020 |
Info@comfortmedical.se | www.comfortmedical.se