

CM ArmFix

CM ProLine



Article no: CMPARM

PRODUCT DESCRIPTION

CM ArmFix is a stable high positioning cushion for forearms. The core of the cushion is made of pressure relieving polyether foam intended for good comfort. The cover is made of durable PU-cover. Two fixation straps are included every cushion.

- **STABLE**
- **EASY TO USE**
- **COMFORTABLE**

Article no: CMPARM

Size: 45x20 cm, Height 37 cm

Cover: PU-cover

Core: Polyether foam

CE

2022-04-01 (Edition:3)

Distributor:

COMFORT medical

Manufacturer: Bima Plastteknik AB,
Trålvägen 5, 43274, Träslövsläge, Sweden

Comfort Medical AB | Hjalmar Brantings väg 11G | 26143
Landskrona | Sweden | +46 (0)733-109020 |
Info@comfortmedical.se | www.comfortmedical.se