

MRI CUSHIONS

CM ProLine



Article no: CPMRO

Article no: CMPMRHS



PRODUCT DESCRIPTION

MRI cushions for comfort and stabilizing the head. The set of four cushions offer good support for the head. The cushions are soft, flexible and easy to clean.

● FLEXIBLE USE

● COMFORTABLE

Article no: CPMRO

Size: 14x9x2-7 cm

Cover: PU-cover

Core: Viscoelastic foam

Article no: CMPMRHS

Size: 14x9x2 cm

Cover: PU-cover

Core: Viscoelastic foam

Cleaning instructions:

Disinfect with alcohol and / or soap solution. Remove blood stains as soon as possible with cold water.

Recycling: No disassembly required. Sort as combustible waste.

Safety instructions:

- Clean the headrest according to given instructions after use.
- Examine the product regularly for damages. Avoid contact with sharp objects and warm sources, it can damage the product. Damaged product should not be in use.

CE

2021-03-07 (Edition:1)

Distributor: **COMFORT** medical

Manufacturer: Bima Plastteknik AB

Comfort Medical AB | Hjalmar Brantings väg 11G | 26143
Landskrona | Sweden | +46 (0)733-109020 |
Info@comfortmedical.se | www.comfortmedical.se