

MRI CUSHIONS

CM ProLine



Article no: CMPMRO

Article no: CMPMRHS



PRODUCT DESCRIPTION

MRI cushions for comfort and positioning of the head. The set of four cushions offer good support for the head. The cushions are soft, flexible and easy to clean. The cushions are recommended to be used in pairs.

● FLEXIBLE USE

● COMFORTABLE

Article no: CMPMRO

Size: 14x9x2-7 cm

Cover: PU-cover

Core: Viscoelastic foam

Article no: CMPMRHS

Size: 14x9x2 cm

Cover: PU-cover

Core: Viscoelastic foam

Cleaning instructions:

Disinfect with alcohol and / or soap solution. Remove blood stains as soon as possible with cold water.

Recycling: No disassembly required.
Sort as combustible waste.

2022-04-01 (Edition:2)

Distributor: **COMFORT** medical

Manufacturer: Bima Plastteknik AB, Trålvägen 5,
43274, Träslövsläge, Sweden

Comfort Medical AB | Hjalmar Brantings väg 11G | 26143
Landskrona | Sweden | +46 (0)733-109020 |
Info@comfortmedical.se | www.comfortmedical.se