

Knee Cushions


CM ProLine



Article no: CMPKK
Name: Knee cushion
Size: 45x31cm
Height: 17 cm
Cover: Artificial leather (PU)
Core: Polyether foam



Article no: CMPKKF
Name: Knee cushion flack
Size: 55x38 cm
Height: 16-8 cm
Cover: Artificial leather (PU)
Core: Polyether foam




Article no: CMPKKC
Name: Cylinder knee cushion
Size: 50x15 cm
Cover: Artificial leather (PU)
Core: Polyether foam

Article no: CMPLKM25
Name: Small knee cushion 25
Size: 25x17cm
Height: 8 cm
Cover: PU-cover
Core: Polyether foam

Article no: CMPLKM40
Name: Small knee cushion 40
Size: 40x17 cm
Height: 8 cm
Cover: PU-cover
Core: Polyether foam

Article no: CMPLKH40
Name: Small knee cushion 40 (Hard)
Size: 40x17 cm
Height: 8 cm
Cover: PU-cover
Core: Sensaten



● **COMFORTABLE**

● **ERGONOMIC**

PRODUCT DESCRIPTION

CM Proline offer a variety of comfortable knee cushions. Main objectives and purpose are hygiene, ergonomics and comfort. Knee cushions are available in a range of shapes to cover all needs of positioning. All knee cushions are covered in outer material that can easily be cleaned. Knee cushions in artificial leather are delivered with ergonomic handles.

Cleaning instructions: Disinfect with alcohol and / or soap solution. Remove blood stains as soon as possible with cold water.

CE

2021-02-10 (Editon:1)

Distributor:

COMFORT medical

Manufacturer: Bima Plastteknik AB, Trålvägen
5, 43274, Träslövsläge, Sweden

Comfort Medical AB | Hjalmar Brantings väg 11G | 26143
Landskrona | Sweden | +46 (0)733-109020 |
Info@comfortmedical.se | www.comfortmedical.se