


EN

## Buttoned Sleeves



MAVIG

## Sewn on Sleeves



Article-No.: **RASLR**  
**RASLF**

Sleeves for all-round protection (RA631 and RA632)  
Sleeves for front protection (RA640, RA641, and RA660)

### Please note:

*The protective surface of the front (surface without pinholes) ends at the seams. This also applies to the sewn-on sleeves.*

*With the button-on sleeves, a continuous protective surface is ensured by the overlap.*

*Made in Germany*

## X-Ray Protective Clothing Optional with Sleeve(s)

For all model ranges (except RA650), you can optionally choose sleeves in order to extend the radiation protection\*. The radiation on the shoulder joint is significantly reduced by the sleeves. Likewise, incident radiation on the side via the armholes is considerably reduced.

The standard lead value is 0.25 mmPb. This ensures good mobility during the daily working routine. The shoulder protector is always manufactured in the colour "Titan".

The correct size as well as the protective material depend on the chosen X-ray protective clothing. MAVIG offers the sleeves both as a button-on solution and firmly sewn-on. You can also choose whether you want only one sleeve, on the right or left, or pairs.

- **Protects the complete AC shoulder joint with standard: 0.25 mmPb**  
(on request: 0.35 mmPb or 0.50 mmPb)
- **Protective material options: NovaLite, Leadfree, or Standard Lead**
- **Outer material in the colour Titan, optional: ComforTex® HPMF or ComforTex® HPMF Hybrid**
- **Connected by snap fasteners or sewed on**
- **Available individually or per pair**
- **Sizes ranging from XS to 3XL (matching the PPE)**

### MAVIG GmbH

X-Ray Protection and  
Medical Suspension Systems

Stahlgruberring 5  
81829 Munich  
Germany

PO Box 82 03 62  
81803 Munich  
Germany

Phone +49 (0) 89 / 420 96 - 0  
Fax +49 (0) 89 / 420 96 - 200  
e-Mail [info@mavig.com](mailto:info@mavig.com)