



<b>Models:</b>	<b>RP2761E</b>	Testicle Protection for men
	<b>RP2761K</b>	Testicle Protection for boys
	<b>RP2761G6232</b>	Testicle Shield Set for Children aged 1 to 10

Made in Germany



## RP276 – Testicle Shield for Men/Boys/Children

This male gonad protection model covers the penis and scrotum all around with a lead equivalent of 1.00 mmPb. Of course, MAVIG has made sure that there is no risk of injury to the patient when using the protection.

The hook-and-loop fastener provides a good hold when adjusted correctly. For cleaning and disinfection, this model series can be unfolded using the pressstuds. We also have suitable disposable hygiene bags\* in our range.

• <b>Protective Material</b>	Standard lead
• <b>Lead Equivalent</b>	1.00 mmPb
• <b>Outer Cover/Colour</b>	Soft PVC cover in light grey (boys/men); Textile cover made of special, easy to clean material in Teddy design (infants/baby/children)
• <b>Sizes</b>	RP2761E - for men (dimensions approx. 155 x 120 x 90 mm) RP2761K - for boys (dimensions approx. 110 x 90 x 65 mm) RP2761G6232 - Set with 3 protections from infant to child sizes

**Hint:** Different hygienic cover sets are available for the sizes E and K, which each contain 100 covers. The set RP-HYG-CEK contains 100 covers for size E and 100 covers for size K.

**MAVIG GmbH**  
X-Ray Protection and  
Medical Suspension Systems

Stahlgruberring 5  
81829 Munich  
Germany

PO Box 82 03 62  
81803 Munich  
Germany

Phone +49 (0) 89 / 420 96 - 0  
Fax +49 (0) 89 / 420 96 - 200  
e-Mail [info@mavig.com](mailto:info@mavig.com)