



Management
System
EN ISO
13485:2016
www.tuv.com
ID 0000075020

Regulation
(EU) 2017/745




Skot Silver

STRETCHER



■ HI-TECH STRETCHER FOR A&E AND PATIENT TRANSFER WITH EXCELLENT INFECTION CONTROL ■

Fully Washable System



Foldable IV pole
fixed to the frame
(optional)



Antibacterial protective casing
with oxygen bottle holder

Brake pedals at both
head and foot sides

Bilateral adjustable height pedals with 'long strokes': less effort is needed to raise the lying surface.



ERGONOMIC

PU cover mattresses for enhanced comfort



PATIENT COMFORT

Fully Washable System: Quickly removable PVC casing and Velcro-free mattress for easy cleaning and enhanced infection control.



HYGIENE

Steering handles and brake pedals accessible from both ends of the stretcher.



HANDLING

Specifications

Overall width (base 620 mm): 740 mm

Overall length: 2100 mm

Adjustable height: 600 to 930 mm

Lying width: 620 / 720 mm

Lying length: 1900 mm

Backrest: 0° to 90°

Trendelenburg/Anti-trend: -15°/+15°

Wheels: ø 200mm with central braking

Mattress thickness: 80 mm

Maximum load: 300 kg



STANDARD FEATURES

Steel frame with white epoxy coating

4 wheels & 5th directional wheel

Central braking

Mattress with PU cover

Standard IV pole

Antichoc bumpers

2 IV pole support

2 foldable barriers

Backrest on gas jack

Fixed handles at head and foot end

Adjustable height pedals on both end

Wheel control pedal at head and foot end

Tray for electric material (optional)



Knee break (optional)

OPTIONAL EXTRAS

(non-exhaustive list
please contact us)

Foldable handles at head end

Foldable handles at foot end

Two support straps

Two stainless steel rails for accessories

Radio upper frame

Radio cassette holder

Knee break

Technical tray with support

Dining tray on barriers

Paper roll holder

Urine pouch holder

Child friendly barriers protections

Non-standard mattress colour

Ultra comfort mattress (pressure point relief)

Easy gliss mattress (allows smoother transfer)

Stomatological headrest

Foldable IV pole fixed to the frame

Telescopic IV pole