



Management
System
EN ISO
13485:2016

www.tuv.com
ID: 0000075020

Directive
93/42/CEE | **CE**

Manchester INN Outdoor

STACKABLE TRANSFER CHAIRS



■ MAXIMISE YOUR STORAGE SPACE ■



IV pole



Coin locking
system" (option)



Anti-theft arch
(option)



AVAILABLE
COLOURS
URBAN

acime
Mindful Design

ACIME FRAME - CONCEPTION AND PRODUCTION SITE

1, avenue de l'Europe - 62720 Rinxent - France +33 (0) 3 21 10 25 30
contact@acimeframe.com - www.acimeframe.com

TRANSFER CHAIRS

Reduce floor space by more than 50% (with 5 stackable chairs) .



ERGONOMIC

Soft upholstery and foldable armrests for easier patient access



PATIENT COMFORT

Hand braking system for increased safety. Exceptional balance and stability.



SAFETY

Large steering handle for super smooth manoeuvrability.



HANDLING

Specifications

Overall width : 705 mm
Overall length : 1125 mm

Seating height : 560 mm
Seating width : 440 mm

Front wheels : Ø125mm

Back wheels : Ø200 mm
Upholstery thickness : 50 mm

Maximum load : 200 kg



Large steering handle with 2 IV pole holders

Passive hand brake and documents holder (option)



Front multidirectional wheels for easier steering

STANDARD FEATURES

Large steering handle
Rilsan-coated steel frame
Passive hand brake
Molded upholstery in polyurethane
Fire proof upholstery M1
Protective bumpers
Two IV pole holders
Foldable foot rest
Two foldable armrests
Belongings holder under the seat
Two lateral hooks



OPTIONAL EXTRAS

*Non-exhaustive list,
please contact us*

Anti-theft arch
IV pole
Safety belt
Documents holder
Non-standard frame colour (blue, yellow, green and red)
Chained wall attachment for locking system
Coin locking system
Key locking system
Personnal identification plate
Urine pouch holder
Oxygen bottle holder
Orthopedic legrest
Orthopedic legrest support only
Manchester INN storage station

