



Management
System
EN ISO
13485:2016

www.tuv.com
ID 000075020

Directive
93/42/EEC | CE

Skot Paediatric

PAEDIATRIC STRETCHER



■ A 360° TRANSPARENT PROTECTION AROUND THE CHILD ■

Fully Washable System



Transparent
panels



Antibacterial
cover

Protection
bumper

Anti pinch system preventing the child, or the nursing staff, to hurt oneself.



ERGONOMIC

Soft mattress with bi-elastic cover. Possibility to add an extra cushion for enhanced comfort (optional).



PATIENT COMFORT

Fully Washable System: PVC casing, Velcro-free mattress and lying surface, removable without any tools in record time.



HYGIENE

The barriers can be moved sideways, or removed at feet and head ends, to give access to any side of the stretcher.



HANDLING

Specifications

Overall width: 710 mm
Overall length: 2120 mm

Adjustable height: 600 to 930 mm
Lying width: 560 mm
Lying length: 1730 mm

Backrest: 0° to 90°
Trendelenburg/Anti-Trend: 15°/15°

Wheels: ø200 mm with central braking

Upholstery thickness: 80 mm

Maximum load: 280 kg

Designed for children care from 3 to 16 years old.



Adjustable backrest up to 90°

STANDARD FEATURES

Steel frame with epoxy coating
4 wheels & 5th directional wheel
Central braking
Mattress with PU cover
Standard IV pole
Protection bumpers
Four IV pole holders
Backrest on gas jacks
Fixed handles at head and feet ends
Bilateral adjustable height pedals
Paediatric barriers with transparent panels
Protection at the head side
Wheel control pedal at the head and foot sides

OPTIONAL EXTRAS

(non-exhaustive list, please contact us)

Paper roll holder
Transparent protection at the head side
Telescopic IV pole
Ultra comfort mattress (memory foam)
Non-standard colour (mattress)



The panels open to provide a complete access to the patient.

IV Pole holder

