



Management
System
EN ISO
13485:2016
www.tuv.com
ID 0000075020


Skot

ACCIDENT & EMERGENCY, SURGERY, PATIENT TRANSFER

Regulation
(EU) 2017/745 | **CE**




■ DURABLE STRETCHER IDEAL FOR SAFE AND SMOOTH PATIENT TRANSFER AND SURGERY ■



Dining tray
(optional)

Fully Washable System



5th wheel and central
braking pedal

Bilateral adjustable
height pedals

acime
Mindful Design

ACIME FRAME - CONCEPTION AND FABRICATION SITE
1, Avenue de l'Europe - 62720 Rinxent - France +33 (0) 3 21 10 25 30
contact@acimeframe.com - www.acimeframe.com/EN/

STRETCHERS

Bilateral adjustable height pedals with 'long strokes': less effort is needed to raise the lying surface.



ERGONOMIC

Mattress with PU cover for enhanced comfort.



PATIENT COMFORT

Fully Washable System: PVC casing, Velcro-free mattress and lying surface, removable without any tools in record time.



HYGIENE

All wheels, including 5th wheel, have high adhesion rubber for enhanced steering.



HANDLING

Specifications

Overall width (base 600mm): 740 mm

Overall length: 2100 mm

Adjustable height: 600 to 930 mm

Available width: 620/720 mm

Lying length: 1900 mm

Backrest: 0° to 90°

Trendelenburg / Anti-trend: -15°/+15°

Wheels: ø 200 mm with central braking

Mattress thickness: 80 mm

Maximum load: 280 kg



Barriers lock out of patient's reach

Barriers with 3 positions

STANDARD FEATURES

Steel frame with epoxy coating
4 wheels & 5th directional wheel
Central braking
Mattress with PU cover
Standard IV pole
Protection bumpers
2 IV pole holders
2 foldable barriers
Backrest on gas jacks
Fixed handles at the foot end
Bilateral adjustable height pedals
Wheel control pedal at the foot end



Backrest on gas jacks

5th wheel for enhance steering

OPTIONAL EXTRAS

(non-exhaustive list, please contact us)

Additional fixed handles at the head end
2 foldable handles at the head end
2 foldable handles at the foot end
Radio-transparent upper frame
X-Ray cassette holder
Wheel control pedal at the head end
Dining tray on barriers
Paper roll holder
Urine pouch holder
"Child-friendly" protection barriers
Stainless universal rail for accessories
Oxygen bottle holder under the backrest
Non-standard colour for the barrier supports
Non-standard colour for the mattress cover
Ultra comfort mattress (pressure relief)
Easy gliss mattress
Stomatology headrest
Articulated IV pole fixed on the frame
Télescopique IV pole

