



Management  
System  
EN ISO  
13485:2016

www.tuv.com  
ID 000075020

# Skiff 3


STRETCHER



Regulation  
(EU) 2017/745 | **CE**



■ LIGHT, EASY TO STEER, ALL CONTROLS ON ONE SIDE FOR HANDLING OPTIMISATION ■




Dining tray  
(optional)



Oxygen bottle  
holder (optional)



Foldable handles  
(optional)



Child friendly barriers  
protection (optional)



Non-standard barrier  
colours (optional)

**acime**  
Mindful Design

ACIME FRAME - CONCEPTION AND PRODUCTION SITE  
1, Avenue de l'Europe - 62720 Rinxent - France +33 (0) 3 21 10 25 30  
contact@acimeframe.com - www.acimeframe.com/EN/

STRETCHERS

Barriers go under the lying surface.  
The mechanism is out of the patient range.



ERGONOMIC

PU cover mattress for enhanced comfort.



PATIENT COMFORT

Velcro-free removable mattress, lying surface and PVC cover.



HYGIENE

All wheels, including 5<sup>th</sup> wheel, have high grip rubber for enhanced steering.



HANDLING

## Specifications

Overall width (620mm): 740 mm  
Overall length: 2065 mm

Adjustable height: 620 to 980 mm  
Lying width: 520/620/720/820 mm  
Lying length: 1900 mm

Backrest: 0° to 90°

Wheels: ø200 mm with central braking

Upholstery thickness : 80 mm

**Maximum load: 215 kg**

Trendelenburg (optional)



Barriers control

## STANDARD FEATURES

Steel frame with epoxy coating  
4 wheels & 5th directional wheel  
Central braking  
Adjustable height pedal at foot end  
Mattress with PU cover  
Standard IV pole  
Protection bumpers  
2 IV pole holders  
2 foldables barriers  
2 fixed handles at foot end  
Backrest on racks

Telescopic IV pole fixed on the frame (optional)



5<sup>th</sup> directional wheel

## OPTIONAL EXTRAS

(non-exhaustive list,  
please contact us)

ø 200 mm outdoor wheels  
Safety straps holder  
Foldable handles at foot end  
Foldable handles at head end  
PU armrest  
Radiotransparent lying surface  
Non-standard colour (mattress)  
Non-standard colour (barriers)  
Fixed/foldable tray for electronic equipment at foot end  
Dinning tray 600x400mm  
Paper roll holder  
Urine pouch holder  
Trendelenburg position  
Child friendly barrier protection  
Stainless steel rail for accessories  
Backrest on gas jack  
Oxygen bottle holder (lower frame)  
Ultra comfort mattress (memory foam)  
Easy gliss mattress (with handles)  
Stomato headrest  
Telescopic IV pole  
Articulated IV pole fixed on frame