





POSITIONING SYSTEMS FOR THE OPERATING ROOM

Bi-Pad

Disposable Anti-Slip Mats
for Trendelenburg Operations



www.biomatrix.it



Bi-Pad
ONEFIT

The classic for optimal
positioning

+NEW



Bi-Pad
DUOFIT

New in bi-coupled
material



+NEW

Bi-Pad DUOFIT

SO02825




EXTRA
GRIP

EXTRA
COMFORT

New technology
two materials
combined together


SO02825 is the latest arrival in BioMatrix family.
A **disposable non-slip pad for the operating room**,
perfect for operations in Trendelenburg position.
This mattress is part of the innovative **BiPad line**.

POSITIONING SYSTEMS FOR THE OPERATING ROOM



BiPad SO02825

composed of two elements
combined together




ULTRA CONFORMABLE

Bi-Pad allows the right combination of pressure, patient position and mobility, in fact, it:

- reduces the risk of injury to the brachial plexus,
- protects the ulnar nerve allowing easy movement of the muscles in the post-operative phase

Disposable anti-slip mat for Trendelenburg operations




Without Bi-Pad
High Pressure
Points



With Bi-Pad
Pressure
reduction



With Bi-Pad DuoFit
reduction and
homogeneous
pressure distribution



+ NEW

Bi-Pad DUOFIT

SO02825

POSITIONING SYSTEMS FOR THE OPERATING ROOM



Medical Device consisting of two distinct layers, designed to **to ensure comfort and safety during surgical operations.**



EXTRA
GRIP

EXTRA
COMFORT

Bi-Pad DUOFIT

SO02825



is made up of two
matched elements



Upper part

VISCOELASTIC MEMORY FOAM

The raw material has production features that ensure excellent grip performance. It perfectly adapts to the patient's body, guaranteeing high adherence and preventing any slipping risks during the operation. It performs an important anti-decubitus function, helping patient's protection from lesions caused by prolonged skin pressure, making the pad highly effective in supporting the patient's well-being during the surgical procedure.



Lower part

POLYURETHANE

The lower part in direct contact with the operating table, made in slow memory PU, guarantees high load-bearing capacity, ensuring homogeneous weight distribution and stable support during surgery. This material, thanks to its gradual adaptation properties, allows minimizing pressure on contact points, contributing to the comfort and safety of the patient during the operations.

Bi-Pad

ONEFIT SO02815




Maximum comfort
for a prolonged
time thanks to its
innovative design

Bi-Pad One Fit is made of **slow-memory viscoelastic polyurethane**, this material allows ergonomic support as well as **complete protection of the patient from injuries caused by skin pressure.**

Bi-Pad has been specifically designed to adapt to the patient's body, its ergonomics allows long-lasting interventions with patients of any age and with a wide range of Body Mass Index

BiPad SO02815

Modular and expandable



BiPad is designed to facilitate anchoring to the operating table during Trendelenburg surgeries, a so-called “anti-shock” position, essential in a wide variety of surgical operations.

Bi-Pad is...



Disposable

It guarantees patient protection in the operating room, eliminating the risk of cross contamination

Customizable

It can be easily combined with other disposable products to facilitate patient positioning

Non-slip

The mat has been studied and designed to facilitate patient's anchoring

Reduction of the costs and times of cleaning and OR preparation

The strength of Bi-Pad is linked to its main production features that allow a reduction in the preparation times of the operating room and an elimination of the costs related to cleaning.

Bi-Pad ONEFIT


SO02815DVG52

BiPad Disposable mat for OR



SO02815S

BiPad disposable mat with silicone for OR



SO02815V

BiPad non-slip disposable mat for OR with Velcro bands and spunlace sheet



SO02815V.I

BiPad non-slip disposable mat for OR with gynecological cut-out, Velcro bands and spunlace sheet



SO02815V.75

BiPad non-slip disposable mat for OR with Velcro bands and spunlace sheet
Size: cm 75x50x3



SO02815VKIT

BiPad non-slip disposable mat for OR with Velcro bands, spunlace sheet and fixing thorax belts




SO02815VKIT.I

BiPad non-slip disposable mat for OR with gynecological cut-out, Velcro bands, spunlace sheet and fixing thorax belts



SO02815VL

BiPad Long non-slip disposable mat for OR with Velcro bands and spunlace sheet
Size: cm 180x52x3



SO02815VKITL

BiPad Long non-slip disposable mat for OR with Velcro bands, spunlace sheet and fixing thorax belts
Size: cm 180x52x3



Bi-Pad DUOFIT


SO02825V

BiPad Coupled disposable non-slip mat for OR with Velcro bands and spunlace sheet



SO02825VKIT

BiPad Coupled disposable non-slip mat for OR with Velcro bands, spunlace sheet and fixing thorax belts



Disposable PU positioners

Code	Size	Dimension	Shape
SO02800	U	ø 23 x 5 cm	
SO02803	U	23x23x6/10	
SO02802	U	23x22x5/10	
SO02805	U	13,5x61x7v	
SO02807	U	42x26,5x14	
SO02808	U	17x18x3,5 cm	
SO02806	U	38x15,5x3,5 cm	
SO02813	U	183x51x5 cm	
SO02811	U	55x35x2,5 cm	
SO02820	U	Ø 15cm - cm 40/60	
SO02823	U	20x23x9.5 cm (± 2%)	
SO02830	U	28x20x17/25	

Supports and bands

Code	Size	Dimension	Shape
SO02819A SO02819B	U	Velcro Band for Legs with Velcro	
SO02814	U	Limb support with Velcro	
SO02815B8	U	Fixing Thorax belt	

BioMatrix

Via Maffini, 23 - 29010 San Giuliano di Castelvetro (PC)
Tel. 0523.826619 Fax. 0523826728
e-mail: info@biomatrix.it

All right reserved. Changes may be made without prior notice.
BioMatrix S.r.l has not to be considered responsible for damages resulting from the lack or inaccuracy of the information reported here.