


Bi-Pad


Disposable Anti-Slip Mat for
Trendelenburg Operations



- + COMFORT
- + MODULARITY
- + FLEXIBILITY

Bi-Pad modular and expandable

Bi-Pad is designed to facilitate **anchoring to the operating table** during operations in Trendelenburg, a so-called "*anti-shock*" position, essential in many surgical operations.




Bi-Pad is made of slow memory viscoelastic polyurethane: this material allows an ergonomic support, protecting patient from injuries caused by skin pressures.




Maximum comfort for a prolonged time thanks to this new design

Bi-Pad has been specially designed to fit patients bodies, its ergonomics allows long-lasting interventions with patients of any age and with a wide Body Mass Index range.

Disposable Anti-Slip Mat for Trendelenburg Operations



Without Bi-Pad
High pressure points




With Bi-Pad
Reduction and distribution of pressure

BI-Pad allows the right combination of pressure, patient position and mobility. Thus,

- it reduces the risk of injury to the branchial plexus,
- it protects the ulnar nerve allowing for easy muscle movement in the post-operative phase.

Bi-Pad is...



Disposable

It guarantees the protection of the patient in the operation room, eliminating the risk of crossed contamination

Customizable

It can be easily combined with other disposable products in order to facilitate patients positioning

Anti-slip

The mat facilitates patients anchoring



Reduction of costs
and both cleaning
and preparation
times

The strength of Bi-Pad is linked to its main features which allow a time and costs reduction in the operating room preparation .


Bi-Pad

At the heart of an integrated positioning system




SOO2815DVG52

- Anti-slip base



SOO2815V

- BI-Pad with straps
- Spunlace tubular sheet




SOO2815V-KIT

- BI-Pad with straps
- Spunlace tubular sheet
- Thorax Belts



SOO2815S

- BI-Pad with silicone



Possibility to create customized kits > Integrated system

Disposable positioning systems

Code	Size	Dimensions	Shape
SOO2800	U	ø 23 x 5 cm	
SOO2803	U	23x23x6/10	
SOO2802	U	23x21x10 cm	
SOO2805	U	13,5x61x7v	
SOO2807	U	42x26,5x14	
SOO2808	U	17x18x3,5 cm	
SOO2806	U	38x15,5x3,5 cm	
SOO2813	U	183x51x5 cm	

The disposable positioners are made of elastic and non-deformable polyurethane foam 25 kgm³ +/-5%, non-toxic and latex free.

Disposable positioning systems

SOO2814	U	Velcro® limb holder	
SOO2814B	U	Limb holder with double Velcro® band	
SOO2815B8	U	Limb holder with double Velcro® band	
SOO2819A	L	Velcro leg bands	
SOO2819B	M		

Holders and bands

BioMatrix

Via Maffini, 23 • 29010 San Giuliano di Castelvetro (PC)

Tel. +39.0523.826619 • Fax 0523.826728 • e-mail: info@biomatrix.it

All rights reserved. Changes may be made without prior notice. BioMatrix S.r.l is not to be considered responsible for damages resulting from the lack or inaccuracy of the information reported here.

